


ข้อมูลด้านเทคนิค

Fluke ES-162P4-2 4-Pole Stake Kit



ภาพรวมผลิตภัณฑ์: Fluke ES-162P4-2 4-Pole Stake Kit

Needed to perform 3- and 4-Pole Fall-of-Potential testing (Stake testing). Used with both the Fluke 1623-2 and Fluke 1625-2. Already included in the Fluke 1625-2 Kit.

รุ่น

**Fluke ES-162P4-2**

- (4) Stakes
- (1) 50 m cable reel of red wire
- (1) 25 m cable reel of blue wire
- (1) 25 m cable reel of green wire

Fluke. ทั่วโลกของคุณ คงอยู่แค่ก้าวต่อไป

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Europe/M-East/Africa
+31 (0)40 267 5100
In Canada (800)-36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com/th-th

©2021 Fluke Corporation. Specifications subject to
change without notice.
11/2021

Modification of this document is not permitted without
written permission from Fluke Corporation.