

LeadWrap Test Lead Wraps



ภาพรวมผลิตภัณฑ์: LeadWrap Test Lead Wraps

LeadWrap fastener for the ability to neatly secure or store test leads and cables.

1. Eliminate tangled leads in your tool box
2. Keeps test leads organized, ready for use
3. Hook and loop fasteners for quick, easy use
4. Three 8-inch (20 cm) test lead wraps per package
5. One year warranty

รุ่น



Fluke LeadWrap

Fluke LeadWrap Test Lead Wrap
Test Lead Wraps, 3 Pk

Fluke. ให้โลกของคุณ คงอยู่แค่ก้าวต่อไป

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Europe/M-East/Africa
+31 (0)40 267 5100
In Canada (800)-36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com/th-th

©2021 Fluke Corporation. Specifications subject to
change without notice.
11/2021

Modification of this document is not permitted without
written permission from Fluke Corporation.