

ข้อมูลด้านเทคนิค

Fluke TL221 SureGrip™ Test Lead Extension Set



คุณลักษณะสำคัญ

- Extends test leads by 1.5 meters
- Silicone-insulated wire resists heat & cold
- Accepts wide range of Fluke test probes and clips

ภาพรวมผลิตภัณฑ์: Fluke TL221 SureGrip™ Test Lead Extension Set

SureGrip™ accessories are designed to improve steadiness in slippery hands. Rubber overmolded surfaces and finger-hugging curves give the user a comfortable, reliable grip on the accessory so they can focus on making an accurate measurement.

- One pair (red, black) of silicone insulated leads with straight connectors on both ends
- Reinforced strain relief
- Includes one pair (red, black) of female couplers
- Extends test leads by 1.5 m
- 600 V CAT IV, 1000 V CAT III, 10 A rating. UL listed
- One year warranty

รุ่น



Fluke TL221

- One pair (red, black) of silicone insulated leads with straight connectors on both ends
- Reinforced strain relief
- Includes one pair (red, black) of female couplers
- Extends test leads by 1.5 m
- 600 V CAT IV, 1000 V CAT III, 10 A rating. UL listed
- One year warranty

Fluke. ให้โลกของคุณ คงอยู่แค่ ละก้าวต่อไป

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Europe/M-East/Africa
+31 (0)40 267 5100
In Canada (800)-36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com/th-th

©2021 Fluke Corporation. Specifications subject to
change without notice.
12/2021

Modification of this document is not permitted without
written permission from Fluke Corporation.