

ข้อมูลด้านเทคนิค

Fluke TL970 Hook and Pincer Set



คุณลักษณะสำคัญ

Handy all-in-one kit that includes 3 Test Lead Sets:

- TL940 Mini-Hook Test Lead Set
- TL950 Mini-Pincer Test Lead Set
- TL960 Micro-Hook Test Lead Set
- One year warranty

รุ่น

**Fluke TL970**

Handy all-in-one kit that includes:

- 3 Test Lead Sets:
 - TL940 Mini-Hook Test Lead Set
 - TL950 Mini-Pincer Test Lead Set
 - TL960 Micro-Hook Test Lead Set

Fluke. ทั่วโลกของคุณ คงอยู่แค่ ปลายนิ้ว

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Europe/M-East/Africa
+31 (0)40 267 5100
In Canada (800)-36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com/th-th

©2021 Fluke Corporation. Specifications subject to
change without notice.
12/2021

Modification of this document is not permitted without
written permission from Fluke Corporation.