

ข้อมูลด้านเทคนิค

## Fluke TLK-220 US SureGrip™ Industrial Test Lead



### คุณลักษณะสำคัญ

Handy all-in-one kit that includes AC220 SureGrip Alligator Clip Set, AC285 SureGrip Large Jaw Alligator Clip Set, TP220 Sharp Test Probes, TL224 SureGrip Silicone Light to Straight Test Leads, Zippered Vinyl Carry Case with Moveable Divider.

### ภาพรวมผลิตภัณฑ์: Fluke TLK-220 US SureGrip™ Industrial Test Lead

Handy all-in-one kit that includes :

- AC220 SureGrip™ Alligator Clip Set
- AC285 SureGrip™ Large Jaw Alligator Clip Set
- TP74 Lantern Tip Test Probes
- TL224 SureGrip™ Silicone Test Lead Set
- Large zippered soft case with moveable divider

รุ่น



## Fluke TLK 220

US SureGrip™ Industrial Test Lead

Included with set:

- AC220 Plunger Style Alligator Clips
- AC285 Large Jaw Alligator Clips
- TP220 Sharp Test Probes
- TL224 Right to Straight Test Leads
- Zippered Vinyl Carry Case with Moveable Divider

**Fluke. ให้โลกของคุณ คงอยู่แค่ ละก้าวต่อไป**

Fluke Corporation  
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:  
In the U.S.A. (800) 443-5853  
In Europe/M-East/Africa  
+31 (0)40 267 5100  
In Canada (800)-36-FLUKE  
From other countries +1 (425) 446-5500  
[www.fluke.com/th-th](http://www.fluke.com/th-th)

©2021 Fluke Corporation. Specifications subject to  
change without notice.  
12/2021

Modification of this document is not permitted without  
written permission from Fluke Corporation.