

ข้อมูลด้านเทคนิค

Fluke TL910 Electronic Test Probes



คุณลักษณะสำคัญ

- Very small 1 mm tips to access hard to reach test points
- Probe tip length up to 100 mm, test lead length: 1 m
- 3 sets of spring loaded gold tips and 2 sets of stainless steel tips included
- CAT II 1000 V, 3 A rating
- One year warranty

ภาพรวมผลิตภัณฑ์: Fluke TL910 Electronic Test Probes

- Very small 1 mm tips to access hard to reach test points
- Probe tip length up to 100 mm, test lead length: 1 m
- Included : 3 sets of spring loaded gold tips and 2 sets of stainless steel tips
- CAT II 1000 V, 3 A rating
- One year warranty

รุ่น



Fluke TL910

Fluke. ให้โลกของคุณ คงอยู่แต่ละก้าวต่อไป

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Europe/M-East/Africa
+31 (0)40 267 5100
In Canada (800)-36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com/th-th

©2021 Fluke Corporation. Specifications subject to
change without notice.
12/2021

Modification of this document is not permitted without
written permission from Fluke Corporation.