

Teknik veriler

TL960 Mikro-Kancalı Test İletkeni Seti



Temel özellikler

- 1 çift (kırmızı, siyah) sarılabilir 4 mm banana fişli ve mikro kancalı test iletkenleri
- 1 mm çapa kadar komponent test iletkenlerine takılabilir mikro kancalar
- 90 cm uzunluğunda PVC yalıtımlı test iletkenleri
- 30 V rms veya 60 V DC, 3 A
- Bir yıl garanti

Ürün İncelemesi: TL960 Mikro-Kancalı Test İletkeni Seti

- 1 çift (kırmızı, siyah) sarılabilir 4 mm banana fişli ve mikro kancalı test iletkenleri
- 1 mm çapa kadar komponent test iletkenlerine takılabilir mikro kancalar
- 90 cm uzunluğunda PVC yalıtımlı test iletkenleri
- 30 V rms veya 60 V DC, 3 A
- Bir yıl garanti

Modeller



Fluke TL960

Micro-Hook Test Leads

Fluke. *Keeping your world up and running.®*

Fluke TURKIYE

P.O. Box 1186
5602 BD Eindhoven
The Netherlands
www.fluke.com.tr

©2022 Fluke Corporation. All rights reserved.

Data subject to alteration without notice.

01/2022

Mer information:

U.S.A. (800) 443-5853 TR +31 (0)40 267 5100 Canada
(905) 890-7600 Övrigt +1 (425) 446-5500

**Modification of this document is not permitted
without written permission from Fluke Corporation.**