

Teknik veriler

Fluke T5-1000 / PRV240 Proving Unit Kit



Modeller



FLUKE-T5-1000-PRV240

1000V Voltage, Continuity and Current Tester

- 4 mm detachable probes
- Two AA batteries (installed)
- Instruction sheet

Fluke PRV240 Proving Unit

- TPAK hanging strap
- Four AA batteries
- Instruction sheet

Fluke. *Keeping your world up and running.®*

Fluke TURKIYE
P.O. Box 1186
5602 BD Eindhoven
The Netherlands
www.fluke.com.tr

©2021 Fluke Corporation. All rights reserved.
Data subject to alteration without notice.
11/2021

Mer information:

U.S.A. (800) 443-5853 TR +31 (0)40 267 5100 Canada
(905) 890-7600 Övrigt +1 (425) 446-5500

**Modification of this document is not permitted
without written permission from Fluke Corporation.**