

Thông số kỹ thuật

Fluke 700PRV Pressure Relief Valve Kit



Tính năng chính

Repeatability $\pm 10\%$ of nominal setting. Multi-turn adjustment screw to set preload on internal disc springs.

Included:

- The 700PRV contains two relief valves
 - One adjustable from 725-2900 psi/50-200 bar (preset at 1360 psi/94 bar)
 - One adjustable from 2900-5800 psi/200-400 bar (preset at 5450 psi/376 bar)
- Instruction sheet

Tổng quan sản phẩm: Fluke 700PRV Pressure Relief Valve Kit

The Fluke 700PRV consists of two relief valves (1360 and 5450 psi) to be used with the 700HTP Hydraulic Test Pump. These relief valves will protect the most commonly used Fluke pressure modules from damage due to accidental over-pressurization. 1/4 BSP male parallel thread to fit Fluke 700HTP.

Thông số kỹ thuật: Fluke 700PRV Pressure Relief Valve Kit

Mechanical and General Specifications

Size	40 x 22 mm
------	------------

Weight	60 g per relief valve
Warranty	One year
Material	Type 303 stainless steel

Model



Fluke 700PRV

Fluke 700PRV Pressure Relief Valve Kit

Fluke. *Giữ cho thế giới của bạn. không ngừng vận động.*

Fluke Corporation

PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:

In the U.S.A. (800) 443-5853

In Europe/M-East/Africa

+31 (0)40 267 5100

In Canada (905) 890-7600

From other countries +1 (425) 446-5500

Representative office of Fluke South East Asia Pte Ltd

C/O Danaher Vietnam

Green Power Tower, 11th Floor Unit 2

35 Ton Duch Thang Street, District 1

Ho Chi Minh City

Vietnam

Tel: +84-8-2220-5371 (ext 103)

Email: info.asean@fluke.com

www.fluke.com/vn

©2022 Fluke Corporation. Specifications subject to change without notice.

01/2022

Modification of this document is not permitted without written permission from Fluke Corporation.