


Thông số kỹ thuật

## Fluke ES-162P3-2 3-Pole Stake Kit



### Tổng quan sản phẩm: Fluke ES-162P3-2 3-Pole Stake Kit

Needed to perform 3-Pole Fall-of-Potential testing (Stake testing). Used with both the Fluke 1623-2 and Fluke 1625-2. This is a subset of what is included in the Fluke 1625-2 Kit.

## Model



### Fluke ES 162P3-2

- (3) Stakes
- (1) 50 m cable reel of red wire
- (1) 25 m cable reel of blue wire

**Fluke.** *Giữ cho thế giới của bạn. không ngừng vận động.*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**  
In the U.S.A. (800) 443-5853  
In Europe/M-East/Africa  
+31 (0)40 267 5100  
In Canada (905) 890-7600  
From other countries +1 (425) 446-5500

**Representative office of Fluke South East Asia Pte Ltd**

C/O Danaher Vietnam  
Green Power Tower, 11th Floor Unit 2  
35 Ton Duch Thang Street, District 1  
Ho Chi Minh City  
Vietnam  
Tel: +84-8-2220-5371 (ext 103)  
Email: info.asean@fluke.com  
www.fluke.com/vn

©2022 Fluke Corporation. Specifications subject to change without notice.  
01/2022

**Modification of this document is not permitted without written permission from Fluke Corporation.**