

Thông số kỹ thuật

Fluke L211 Probe Light Kit

Tổng quan sản phẩm: Fluke L211 Probe Light Kit

Special value combo pairs the new L200 Probe Light with one of the most popular test leads by Fluke and a handy carrying case.

Model

Fluke L211

Fluke L211 Probe Light Kit

L211 Probe Light Kit includes:

- L200 Probe Light, TL71 Premium DMM Test Lead Set and C75 Case
 - Small, rugged light easily attaches to any Fluke test probe
 - Flexible silicone test leads are heat and cold resistant
 - Zippered case has two inside pockets
 - Probe light frees both hands for work
 - 120 hours of battery life
 - Uses 2 3-volt watch batteries
-

Fluke. *Giữ cho thế giới của bạn. không ngừng vận động.*

Fluke Corporation

PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:

In the U.S.A. (800) 443-5853

In Europe/M-East/Africa

+31 (0)40 267 5100

In Canada (905) 890-7600

From other countries +1 (425) 446-5500

Representative office of Fluke South East Asia Pte Ltd

C/O Danaher Vietnam

Green Power Tower, 11th Floor Unit 2

35 Ton Duch Thang Street, District 1

Ho Chi Minh City

Vietnam

Tel: +84-8-2220-5371 (ext 103)

Email: info.asean@fluke.com

www.fluke.com/vn

©2022 Fluke Corporation. Specifications subject to change without notice.

01/2022

Modification of this document is not permitted without written permission from Fluke Corporation.