

Thông số kỹ thuật

VPS101 – Voltage Probe Set, 1:1, one black



Tính năng chính

- 300 V CAT III rated
 - 1:1 Voltage Probe
 - Hook Clip for Probe Tip
 - Ground Lead with Mini Alligator
 - Insulation Sleeve
- 30 V CAT II rated Ground Spring for Probe Tip

Tổng quan sản phẩm: VPS101 – Voltage Probe Set, 1:1, one black

Double insulated voltage probe for ScopeMeter 190 series.

Thông số kỹ thuật: VPS101 – Voltage Probe Set, 1:1, one black

Specification	
Attenuation	1x
Bandwidth	30 MHz
Max. input voltage	300 V CAT III
Max. voltage to ground	600 V CAT III per IEC / EN61010-031

Input resistance	1 M Ω
Input capacitance	42 pF + input capacitance of instrument
Cable length	1.2 m
Warranty	One year

Fluke. *Giữ cho thế giới của bạn. không ngừng vận động.*

Fluke Corporation

PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:

In the U.S.A. (800) 443-5853

In Europe/M-East/Africa

+31 (0)40 267 5100

In Canada (905) 890-7600

From other countries +1 (425) 446-5500

Representative office of Fluke South East Asia Pte Ltd

C/O Danaher Vietnam

Green Power Tower, 11th Floor Unit 2

35 Ton Duch Thang Street, District 1

Ho Chi Minh City

Vietnam

Tel: +84-8-2220-5371 (ext 103)

Email: info.asean@fluke.com

www.fluke.com/vn

©2022 Fluke Corporation. Specifications subject to change without notice.

01/2022

Modification of this document is not permitted without written permission from Fluke Corporation.