

Thông số kỹ thuật

P5544 Two Gauge Stand

Tổng quan sản phẩm: P5544 Two Gauge Stand

This adapter mounts directly to the test port of the calibrator and allows for the calibration of two instruments at the same time or the connection of a reference test instrument. Maximum working pressure is 10,000 psi (700 bar).

Thông số kỹ thuật: P5544 Two Gauge Stand

Specifications	
Body Material	Brass and Stainless steel
O-Rings	Nitrile
Test Port Adapters	1 set, 1/8, 1/4, 3/8 and 1/2 NPT or BSP and adapters supplied with calibrator
Calibrator Test Connection	1/2 BSP male (mating adapter supplied)

Model

P5544

P5544 Two Gauge Stand

Two Gauge Stand

Fluke. *Giữ cho thế giới của bạn. không ngừng vận động.*

Fluke Corporation

PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:

In the U.S.A. (800) 443-5853

In Europe/M-East/Africa

+31 (0)40 267 5100

In Canada (905) 890-7600

From other countries +1 (425) 446-5500

Representative office of Fluke South East Asia Pte Ltd

C/O Danaher Vietnam

Green Power Tower, 11th Floor Unit 2

35 Ton Duch Thang Street, District 1

Ho Chi Minh City

Vietnam

Tel: +84-8-2220-5371 (ext 103)

Email: info.asean@fluke.com

www.fluke.com/vn

©2022 Fluke Corporation. Specifications subject to change without notice.

01/2022

Modification of this document is not permitted without written permission from Fluke Corporation.